

Crimp head and bottle tool for FS10 and FS32




Sales manual

Rev.:03

25.02.2011

AEF



Part number	Description	Material
89-100-051	Bottle tool for FS10, Ø14 mm – 38 mm	POM
89-100-052	Bottle tool for FS10, Ø35 mm – 63 mm	POM
89-100-053	Bottle tool for FS10, Ø57 mm – 73 mm	POM
89-100-050	Bottle tool - special size for FS10	POM
89-100-105	Crimp head DIN 13 mm standard cap	-
89-100-106	Crimp head DIN 13 mm flip-off	-
89-100-107	Crimp head DIN 20 mm standard cap	-
89-100-108	Crimp head DIN 20 mm flip-off	-
89-100-109	Crimp head DIN 32 mm std/flip-off (only for FS32)	-
89-100-110	Crimp head for special size of crimp cap	-
89-120-051	Bottle tool for FS32, Ø20 mm – 51 mm	POM
89-120-052	Bottle tool for FS32, Ø45 mm – 76 mm	POM
89-120-053	Bottle tool for FS32, Ø70 mm – 97 mm	POM
89-120-050	Bottle tool - special size for FS32	POM
89-100-055	Pocket for vials with height less than 55mm or a vials with a very small diameter	POM

- For each size and type of crimp cap an individual crimping head is applied. It is made to suit the specific cap.
- Three standard bottle tools are used for bottle diameters as described above. For other shapes of bottles, a bottle tool is specially made to suit the bottle in question.
- With the order, please send samples of vials stoppers and caps to Flexicon for testing prior to shipment.